

Biochemicals and Reagents for Life Science Research

ALPHABETICAL
LIST

BIOACTIVE
PEPTIDES

IMMUNO-
CHEMICALS

MOLECULAR
BIOLOGY

RBI,
NEUROSCIENCE,
SIGNAL
TRANSDUCTION

TISSUE
CULTURE

OTHER
PRODUCT
GROUPS/USP

EQUIPMENT,
BOOKS AND
SUPPLIES

DIAGNOSTIC
KITS AND
REAGENTS

PRODUCT
INDEX



SIGMA®

1999

ALPHABETICAL
LIST

BIOACTIVE
PEPTIDES

IMMUNO-
CHEMICALS

MOLECULAR
BIOLOGY

RBI,
NEUROSCIENCE,
SIGNAL
TRANSDUCTION

TISSUE
CULTURE

OTHER
PRODUCT
GROUPS/USP

EQUIPMENT,
BOOKS AND
SUPPLIES

DIAGNOSTIC
KITS AND
REAGENTS

PRODUCT
INDEX



P.O. Box 14508, St. Louis, MO 63178 USA

ORDER / CUSTOMER SERVICE 800-325-3010

To place an order or to inquire about an existing order.

Fax: 800-325-5052 24 hours a day, seven days a week.

For large scale inquiries, see page 16.

TECHNICAL SERVICE 800-325-5832

To discuss applications or product-specific questions.

GENERAL INFORMATION 800-521-8956

**NEW
Super
Site** www.sigma-aldrich.com

The SIGMA-ALDRICH Family

 **SIGMA**

 **ALDRICH**


 **Fluka**

 **Riedel-de Haën**

 **SUPELCO**

ALPHABETICAL LIST OF COMPOUNDS

etic acid; N-Methylglycine)



nt: <0.1%
water, 20°C):
55
<0.1%

Mg: <0.0005%
Na: <0.05%
NH₄⁺: <0.05%
P: <0.0005%
Pb: <0.001%
Zn: <0.0005%

NO₂ FW 89.09

NO₂ FW 89.09
10 g 7.65
25 g 13.75
100 g 21.70

NO₂ • HCl
10 g 12.45
25 g 19.45
100 g 71.35

DRIDE
3-dioxopiperazine) 10 g 47.20
25 g 103.80

yellow crystals.
H₁₀N₂O₂ FW 142.2

ROGENASE 10 units 66.10
[receptor]
[demethylating]; EC 1.5.99.1)
[on] as sp.
[under] containing approx. 30%
[balance] approx. 60% sucrose,
[phosphate] buffer salts and trace

units per mg protein.
One unit will convert 1.0 μmole of
[line] and formaldehyde per min at

L ESTER 10 g 31.00

C₂H₁₁NO₂ • HCl FW 153.6

ASE
[gen] oxidoreductase (demethylating);

One unit will form 1.0 μmole of
[m] sarcosine per min at pH 8.3 at

cter species 25 units 22.45
wder 100 units 61.65
rox. 50% 250 units 114.95

[balance]
ssium gluconate and EDTA.
nits per mg protein.

species 1,000 units 91.50
wder
units per mg solid.

acterium 25 units 48.95

wder containing approx. 50%
[balance] primarily phosphate
tose.
nits per mg protein.

(Continued)

PRODUCT
NUMBER

US \$

(Continuation of)
SARCOSINE OXIDASE

S 5896 From *Pseudomonas* species 500 units 105.00
Lyophilized powder
containing approx. 30% protein (Biuret).
Activity: Approx. 5 units per mg protein.

SARCOSYL-

Peptides are listed using standard 3-letter
abbreviations. For N-terminal sarcosyl-peptides see
SAR- beginning on Page 933

SARKOSYL

See: N-Lauroylsarcosine, Sodium Salt Page 626

SAR-PRO-ARG p-NITROANILIDE 5 mg 44.75
S 9009 Dihydrochloride 10 mg 74.30
Chromogenic substrate for throm-
bin. 25 mg 147.75
Ref.: Duncan, A., et al., Clin. Chem., 31, 853 (1985).
[75241-23-5] C₂₀H₃₀N₆O₅ • 2HCl FW 535.4

SARSASAPOGENIN 25 mg 36.80
S 8534 ([25S]-Spirostan-3β-ol; Parigenin) 100 mg 101.95
Minimum 98%
Sarsasapogenin and smilagenin are epimers at the
25-position. See: Page 945
[126-19-2] C₂₇H₄₄O₃ FW 416.6

SASRIN RESINS

See: FMOc-Amino Acid Sasrin Esters Page 457

SATA


See: S-Acetylthioglycolic Acid N-Hydroxysuccinimide
Ester Page 44

SAUVAGINE

See: Bioactive Peptides Page 1162

SAXITOXIN

S 1417 (STX) 10 μg 41.90
Diacetate Salt
From *Protogonyaulax* sp.
Minimum 95% (TLC)
Neurotoxin produced by dinoflagellates in "red tides".
Shellfish consume these dinoflagellates and become
toxic. Mode of action is as sodium channel blocker.
Not assayed by Sigma.
Ref.: 1. Levin, R.E., J. Food Chem., 15, 405 (1991).
2. Narahashi, T., Fed. Proc., 31, 1124 (1972).



[35523-89-8] C₁₀H₁₇N₃O₄ • 2C₂H₄O₂ FW 419.4
R: 26/27/28 S: 45-26-36/37/39-22

neo-SAXITOXIN 10 μg 32.60
S 0170 (neo-STX) Shipped in dry ice
From *Dinoflagellates*
Solution in 0.03 N acetic acid, pH 3.
One of a group of neurotoxins produced by the
organism responsible for "red tides".
Sodium channel blocker
Ref.: 1. Yieytes, M.R., et al., Anal. Biochem., 221, 87
(1993).
2. Wichmann, C.F., et al., Tetrahedron Lett., 22,
1941 (1981).
[64296-20-4] C₁₀H₁₇N₃O₅ FW 315.3
R: 26/27/28 S: 45-26-36/37/39-22

PRODUCT
NUMBER

US \$

NEW SB 203580 1 mg 85.00
S 8307 (4-[4-Fluorophenyl]-
2-[4-methylsulfonylphenyl]-5-[4-pyridyl]-1H-imidazole)
Specific inhibitor of p38 MAP kinase/reactivating
kinase.
Ref.: Ward, S.G., et al., Biochem. Soc. Trans.,
25(2), (1997).

SB(3-n)

This nomenclature refers to alkyl sulfobetaine
detergents, N-alkyl-N,N-dimethylammonio-
1-propanesulfonates. See:
N-Octyl-... (SB3-8) Page 772
N-Decyl-... (SB3-10) Page 327
N-Dodecyl-... (SB3-12) Page 400
N-Tetradecyl-... (SB3-14) Page 990
N-Hexadecyl-... (SB3-16) Page 540
N-Octadecyl-... (SB3-18) Page 771

SBFI

S 0762 Tetraammonium Salt 500 μg 143.85
1 mg 223.50
Approx. 90%
Fluorescent indicator for intracellular sodium
Ref.: 1. Minta, A., and Tsien, R.Y., J. Biol. Chem.,
264, 19449 (1989).
2. Chen, X., and Gross, R.W., Biochem., 33, 13769
(1994).
[124549-08-2] C₂₄H₃₂O₁₃N₆ FW 906.9

SBFI-AM

S 1148 Cell permeant sodium selective 100 μg 44.95
fluorescent indicator 1 mg 249.40
Ref.: 1. Minta, A., and Tsien, R.Y., J. Biol. Chem.,
264, 19449 (1989).
2. Chen, X., and Gross, R.W., Biochem., 33, 13769
(1994).
[129423-53-6] FW 1127.1

SCANDIUM OXIDE

30,787-4 Rare earth content is 99.9% 1 g 73.80
minimum, expressed as Sc₂O₃. 5 g 267.70
Not assayed by Sigma.
Aldrich Brand. Formerly Sigma Product S 1625.
[12060-08-1] Sc₂O₃ FW 137.9

R-(+)-SCH-23390

C 0206 (R-(+)-CHMB: R-(+)-7-Chloro- 5 mg 47.00
8-hydroxy-3-methyl-1-phenyl- 10 mg 75.00
2,3,4,5-tetrahydro-1H-3-benzazepine)
Hydrochloride
Selective D₂ dopamine receptor antagonist
Ref.: O'Boyle, et al., J. Neurochem., 48, 1039
(1987).
[125941-87-9] C₁₇H₁₈ClNO • HCl FW 324.2

SCHAEGLER AGAR

See: Microbiological Media and Components
Page 2138

SCHARDINGER α-DEXTRIN

See: α-Cyclodextrin Page 311

SCHARDINGER β-DEXTRIN

See: β-Cyclodextrin Page 311

SCHENK AND HILDERBRANDT BASAL SALT
MIXTURE

See: Plant Tissue Culture Media and Reagents
Page 1890

SCHENK AND HILDERBRANDT VITAMIN MIXTURE

See: Plant Tissue Culture Media and Reagents
Page 1890

ALPHABETICAL LIST OF COMPOUNDS

PRODUCT NUMBER	US \$	PRODUCT NUMBER	US \$
L 0631 (Behenyl laurate) Approx. 99% [42231-82-3] C ₂₂ H ₄₄ O ₂ FW 508.9	25 mg 8.25	L 5503 (Dodecanoyl chloride) Approx. 99% d = 0.92 g/ml Sealed ampule. [112-16-3] C ₁₂ H ₂₅ ClO FW 218.8 R: 34 S: 26-27-36/37/39-3/7/9	100 ml 25.20
L 4500 (Butyl laurate) Approx. 98% (capillary GC) d = 0.86 g/ml [106-18-3] C ₁₆ H ₃₂ O ₂ FW 256.4	100 ml 16.40	L 7262 (25234-60-0) C ₁₇ H ₃₅ NO ₂ Cl FW 321.9	10 g 33.45
L 2136 (N,N-Dimethylauramide) Approx. 95% d = 0.87 g/ml [3007-53-2] C ₁₄ H ₂₉ NO FW 227.4	100 g 51.00	LAUROYL COENZYME A See: Coenzyme A and Derivatives Page 283	
L 4625 (Ethyl laurate) Minimum 99% d = 0.86 g/ml [106-33-2] C ₁₄ H ₂₈ O ₂ FW 228.4	50 ml 7.70 100 ml 12.15 500 ml 40.40	LAUROYL-(CARBONYL-¹⁴C) COENZYME A See: Radiochemicals Page 283	
L 3900 (N-Lauroyloxysuccinimide) Approx. 98% (TLC) Intermediate in the synthesis of lauroyl derivatives of amines and thiols. Ref.: Al-Arif, A., et al., J. Lipid Res., 10, 344 (1969). Lapidot, Y., et al., Biochim. Biophys. Acta, 145, 292 (1967). [14565-47-0] C ₁₈ H ₃₇ NO ₄ FW 297.4	1 g 13.90 5 g 45.50	1-LAUROYL-2-OLEOYL-3-PAL-MITOYL-rac-GLYCEROL (C12:0/C18:1[cis]-9/C16:0) (1-Dodecanoyl-2-[cis-9-octadecenyl]-3-hexadecanoyl-rac-glycerol) Approx. 98% [51604-51-4] C ₄₈ H ₉₂ O ₆ FW 777.3	25 mg 117.30
L 4750 (Lauryl laurate) [13945-76-1] C ₂₄ H ₄₈ O ₂ FW 368.6	5 g 14.10 100 g 90.30	N-LAUROYL-SARCOSINE $\text{CH}_3(\text{CH}_2)_9\text{CH}_2 - \overset{\text{O}}{\parallel} \text{C} - \text{N}(\text{CH}_3) - \overset{\text{O}}{\parallel} \text{C} - \text{OH}$	
L 7272 (Methyl laurate) 99+% (capillary GC) d = 0.87 g/ml Also available as part of a kit. See: Standards and Controls Section Page 2199 [111-82-0] C ₁₃ H ₂₆ O ₂ FW 214.3	1 g 5.85 5 g 19.35 25 g 61.50	L 5000 Free Acid Minimum 95% May be liquid at room temperature. [9778-9] C ₁₅ H ₃₁ NO ₃ FW 271.4 R: 36/37/38 S: 26-36	50 g 10.90 100 g 16.20 500 g 52.20
L 9631 (Myristyl laurate) Approx. 99% [22412-97-1] C ₂₆ H ₅₂ O ₂ FW 396.7	100 mg 17.20	L 5777 Sodium Salt SigmaUltra Minimum 97% Solubility (0.1 M in water, 20°C): complete, very faint yellow Insoluble matter: <0.1% Al: <0.0005% Ca: <0.005% Cu: <0.0005% Fe: <0.005% K: ≤100 ppm [137-16-6] C ₁₅ H ₃₁ NO ₃ Na FW 293.4	50 g 28.20 100 g 46.50 500 g 156.10
L 9255 (Oleoyl laurate) Approx. 99% Sealed ampule (liquid). [19149-85-0] C ₃₀ H ₅₈ O ₂ FW 450.8	1 g 57.50	L 5125 Sodium Salt Minimum purity: 94% See also: Molecular Biology Products Page 1543 [137-16-6] C ₁₅ H ₃₁ NO ₃ Na FW 293.4	50 g 12.50 100 g 21.30 500 g 66.40 1 kg 118.50
L 8641 (Palmityl laurate) Approx. 99% [20834-06-4] C ₂₈ H ₅₆ O ₂ FW 424.7	100 mg 16.75 1 g 53.05	L 5125 Sodium Salt Minimum purity: 94% See also: Molecular Biology Products Page 1543 [137-16-6] C ₁₅ H ₃₁ NO ₃ Na FW 293.4	50 g 12.50 100 g 21.30 500 g 66.40 1 kg 118.50
L 3253 (Propyl laurate) Approx. 99% d = 0.86 g/ml [3681-78-5] C ₁₅ H ₃₀ O ₂ FW 242.4	5 g 38.45	L 4003 (Lauryl acetate) Minimum 99% d = 0.86 g/ml [112-66-3] C ₁₄ H ₂₈ O ₂ FW 228.4	1 g 19.20
L 1632 (Dodecanoic anhydride) Approx. 99% [645-66-9] C ₂₄ H ₄₆ O ₃ FW 382.6 R: 36/37/38 S: 26-36	1 g 15.70 10 g 81.40 25 g 146.60	L 5375 (1-Dodecanol; Dodecyl alcohol) Approx. 99% Also available as part of a kit. See: Standards and Controls Section Page 2200 [112-53-8] C ₁₂ H ₂₆ O FW 186.3 R: 36/37/38 S: 26-36	500 g 11.90
L 3131 α-LAUROYL-CARNITINE CHLORIDE [7023-03-2] C ₁₉ H ₃₈ NO ₄ Cl FW 380.0	25 mg 14.05 100 mg 36.45 250 mg 71.20 1 g 194.75	LAURYLAMINE See: Dodecylamine Page 400	
		LAURYL BROMIDE See: 1-Bromododecane Page 195	
		1-O-LAURYL-rac-GLYCEROL See: 1-O-Dodecyl-rac-glycerol under Lipids Page 401	

PRODUCTS

PURIFICATION

GER RNA ISOLATION KITS
ard Isolation Kit 1 kit 255.45
 ent for 6 preparations from
 10⁶–3 × 10⁶ cultured cells or 400–1,000 mg
 ie
 22-42-36/37/38-41 S: 26-36

REAGENTS

3-CHLOROPROPANE 200 ml 19.80
 e for use in RNA extractions using any of the
 agents. BCP can be used in place of
 form, and is less toxic than chloroform. It does
 adversely affect quality or quantity of the isolated
 . Chomczynski, P. and Mackey, K., Anal.
 m., 225, 163 (1995).
 [67-63-0] C₃H₇BrCl FW 157.4
 16/37/38 S: 26-36

OL
 er: Molecular Biology Reagents Page 1623

CHLORIDE
 [78] CsCl FW 168.4
 99+% 5 g 11.30
 30% solution <0.1 25 g 38.60
 50% solution <0.02 100 g 98.50
 11% 500 g 390.10
 01% 6 × 100 g 472.70
 RNase and 1 kg 571.60
 e - None detected. 6 × 500 g 1872.00
 >98% 25 g 16.80
 30% solution <0.2 50 g 28.00
 RNase and protease - 100 g 46.60
 etected. 250 g 102.55
 500 g 184.80
 1 kg 276.15

FORM 25 ml 7.65
 99+% 4 × 25 ml 23.45
 for use in nucleic acid 500 ml 25.30
 ion. When used 6 × 500 ml 121.10
 r with phenol, the
 cy of protein extraction from crude DNA is
 d. Can also be used to remove traces of
 from aqueous DNA and RNA samples.
 d with amylenes.
 mbrook, J., et al., Molecular Cloning: A
 ry Manual, Cold Spring Harbor Laboratory
 p. E.3-E.4.
 [77-47-6] CHCl₃ FW 119.4
 3-40-48/20/22 S: 36/37

ORM: ISOAMYL ALCOHOL 1 pt 25.50
 1 qt 41.10
 ed with amylene (2-methyl-
 ie)
 for use in the purification of nucleic acids
 mbrook, J., et al., Molecular Cloning: A
 ry Manual, Cold Spring Harbor Laboratory
 p. E.3-E.4.
 5-23/24/25-41 S: 45-26-36/37/39

3ONUCLEASE I
 ar: Cloning, Modifying Enzymes Page 1556

a-aldrich.com/techinfo

MOLECULAR BIOLOGY PRODUCTS

NUCLEIC ACID EXTRACTION AND PURIFICATION

ENDOTOXIN REMOVAL SOLUTION 25 ml 25.00
E 4274
 A solution of Triton X-114 in Tris-buf-
 fered saline, pH 7.4. Useful for removal of endotoxins
 from a solution of plasmid DNA by phase separation.¹
 A 25 ml aliquot is sufficient to treat 500 ml of lysate,
 from which 8-10 mg of endotoxin-free DNA can be
 isolated.
 Ref.: Cotton, M. et al., Gene Ther., 1, 239 (1994).
 R: 36/37/38 S: 26-36

ETHANOL
 (Ethyl alcohol)
 Suitable for use in the precipitation of nucleic acids.
 US excise tax included; ATF license not required.
 [64-17-5]
 R: 11 S: 7-16

Absolute (200 proof) 500 ml 25.50
E 7023
 Water ≤0.005%
95+% (190 proof) 1 gal 110.00
E 7148

GUANIDINE THIOCYANATE
 See under: Molecular Biology Reagents Page 1625

HEXADECYLTRIMETHYLAMMONIUM BROMIDE
 See under: Molecular Biology Reagents Page 1625

ISOAMYL ALCOHOL 25 ml 7.15
 (Isopentyl alcohol; 3-Methyl- 4 × 25 ml 21.90
 1-butanol) 500 ml 23.60
 ≥98.5% 6 × 500 ml 112.95
 Suitable for use in nucleic acid
 purification.
 Ref.: Sambrook, J., et al., Molecular Cloning: A
 Laboratory Manual, Cold Spring Harbor Laboratory
 (1989) p. E.3-E.4.
 [123-51-3] C₅H₁₂O FW 88.15
 R: 10-22-41/37/38 S: 16-26-36/37/39-23

ISOPROPANOL 25 ml 7.25
 (2-Propanol) 4 × 25 ml 22.05
 Purity: 99+% 500 ml 23.75
 Water ≤0.05%
 Suitable for use in the precipitation of nucleic acids.
 When compared to ethanol, 50% less is required for
 nucleic acid precipitation, thus minimizing the total
 volume to be centrifuged for DNA or RNA recovery.
 Ref.: Sambrook, J., et al., Molecular Cloning: A
 Laboratory Manual, Cold Spring Harbor Laboratory
 (1989) p. E.13-E.14.
 [67-63-0] C₃H₈O FW 60.10
 R: 11 S: 7-16

N-LAUROYLSARCOSINE 50 g 15.35
L 9150
 Sodium Salt 100 g 23.85
 Purity: ≥97% 250 g 47.75
 DNase and RNase: None detected.
 Soluble in concentrated salt solutions that are com-
 monly used in the cell lysis step of RNA purification
 Ref.: Sambrook, J., et al., Molecular Cloning: A
 Laboratory Manual, Cold Spring Harbor Laboratory, p.
 7.20 (1989).
 [137-16-6] C₁₃H₂₇NO₃Na FW 293.4

LITHIUM CHLORIDE
 See under: Molecular Biology Reagents Page 1626

LYSOZYME
 See under: Cloning, Modifying Enzymes Page 1557

PHASE DIVIDER GEL

For organic extractions of nucleic acids—forms a
 durable barrier between phases, facilitating decanting
 or pipetting into a fresh tube.
 • Improve nucleic acid recovery 20-30%
 • Reduce contamination from interphase
 • Gel is inert; does not interfere with subsequent
 enzymatic reactions
 • Convenient for samples from 10 µl to 20 ml
 The gel comes pre-dispensed in common sizes of
 centrifuge tubes. Just add your phenol or
 phenol/chloroform extraction mixture and spin. The
 gel migrates under centrifugal force, forming a
 barrier above the organic phase and interphase.
Usage Suggestions

Organic phase →	phenol: chloroform: isoamyl alcohol	chloroform: isoamyl alcohol	water- or buffer-sat'd phenol: chloroform	water- or buffer-sat'd phenol
↓ Aqueous phase				
<0.5 M NaCl	L,H	L,H	L,H	L
<1 mg/ml BSA	L,H	L,H	L,H	L
Cleared lysate	H	H	H	NC
Tissue homogenate	L,H	L,H	L,H	L
RNA purification	H	H	H	NC

Notes: L = Light; H = Heavy; NC = not compatible

P 1723 Phase Divider Light/1.5 6 / pkg 10.40
 in 1.5 ml microcentrifuge tubes 200 / pkg 109.20
 for 50-500 µl samples

P 2098 Phase Divider Light/15 10 / pkg 12.50
 in 15 ml low-speed centrifuge 50 / pkg 48.90
 tubes for 1-6 ml samples 100 / pkg 87.40

P 2348 Phase Divider Light/50 10 / pkg 32.25
 in 50 ml centrifuge tubes for 5-20 ml samples

P 2598 Phase Divider Light/Syringe 10 / pkg 59.30
 in 3 ml syringes for dispensing into any size centrifuge tube

P 1848 Phase Divider Heavy/1.5 6 / pkg 10.40
 in 1.5 ml microcentrifuge tubes 200 / pkg 109.20
 for 50-500 µl samples

P 1973 Phase Divider Heavy/10 100 / pkg 87.40
 in 10 ml high-speed centrifuge tubes for 1-4 ml samples

P 2223 Phase Divider Heavy/15 10 / pkg 12.50
 in 15 ml low-speed centrifuge 50 / pkg 48.90
 tubes for 1-6 ml samples 100 / pkg 87.40

P 2473 Phase Divider Heavy/50 10 / pkg 32.25
 in 50 ml centrifuge tubes for 5-20 ml samples

P 2723 Phase Divider Heavy/Syringe 10 / pkg 59.30
 in 3 ml syringes for dispensing into any size centrifuge tube

PHENOL 25 g 26.25
P 1037 99+% 100 g 38.85
 Redistilled and packaged under 500 g 88.20
 nitrogen in amber bottles. 6 × 100 g 186.50
 Contains no preservatives
 Mg ≤2 ppm
 Pb ≤5 ppm
 Melting point: 40 ± 2°C
 [108-95-2] C₆H₆O FW 94.11
 R: 24/25-34 S: 28-45

Molecular Biology Products

To place an order call 800-325-3010 • www.sigma-aldrich.com/order

1543